Pseudomonas aeruginosa VNTR typing form						
Effective Date	18/04/2023	Document Type	Request Form			
Active Date	18/04/2023	Page Number	Page 1 of 1			
Author/Reviewer	L Crawford	Authoriser	L McCoy			
Revision Number	3.0	Document Number	M-2057			

Send to: Molecular Typing, Microbiology, Kelvin Building, Royal Hospitals Trust, Grosvenor Road, Belfast BT12 6BA. Direct Tel 028 96 151646 (9am-5pm Mon-Fri). Email: Moleculartyping@belfasttrust.hscni.net

VEELA I	ARFI	OPI	ENITED	DETAILS	I ECIRI V
ALLIAL	ADLL			DL I AILO	LLGIDLI

Male/Female	Surname	Forename	D.O.B.			
H&C No						
Hospital	Ward/Clinic					
Or other source						
Referring Laboratory:						
Referring Laboratory number:						
Date and time of sample: Date of referral:						
Site of Isolate/specimen type:						
URGENT* □ ROUTINE □ Details of isolates linked to Urgent requests:						
*Urgent typing is available for isolates from neonates, or on request from a consultant in other cases. Urgent typing of single isolates is not justified unless there is at least one epidemiologically linked clinical or environmental isolate for comparison – provide details of linked samples if urgent testing is required. IF AN ISOLATE IS REFERRED FOR URGENT TYPING, PLEASE ADVISE US BY TELEPHONE OR EMAIL SO THAT RECEIPT CAN BE CONFIRMED						
Surveillance typing. PHA criteria for surveillance reporting and typing of <i>P. aeruginosa</i> isolates are (i) colonisation OR infection of a neonate, or (ii) isolation from a sterile site in a non-neonatal augmented care setting. Isolates for VNTR typing should be sent as freshly grown subcultures on NA slopes. Note: The VNTR typing scheme is specific to <i>Pseudomonas aeruginosa</i> . Other <i>Pseudomonas</i> spp. (i.e. <i>putida</i> , <i>fluorescens</i> , <i>stutzeri</i>) are non-typeable.						
	& time received:	Lab Use:				

